

Scanning Fabry-Perot Interferometers (High-Resolution Spectrometers)

NEW
versions



SA200

Ultra-Stable
Invar Cavity



SA210

Features

- Ultra-Stable Invar Cavity
- Confocal Cavity Design
- Mounted Input and Output Alignment Irises

Specifications

- 1.5 GHz or 10 GHz Free Spectral Range
- Low Scan Voltage, 5 V per FSR @ 633 nm
- 13 Models for the 370 - 1625 nm Range

The SA200 family includes 13 Fabry-Perot interferometers that cover six spectral regions with either a 1.5 GHz or 10 GHz Free Spectral Range (FSR). The design of the Fabry-Perot interferometer cavity is comprised of an Invar cavity with internal piezo stacks. This design utilizes the negative thermal coefficient of the piezo stacks to create the nearly athermal cavity that is necessary for the stability of these high-resolution spectrometers.

The tutorial on pages 810-811 covers the basic theory of operation, including an explanation of the effect of the input beam diameter on the resolution of the instrument. In order for the instrument to achieve the specified resolution, the input beam diameter must not exceed the maximum diameter specification, even though the input aperture for the instrument is significantly larger than the specification.

The SA200 can be mounted via a Ø2" ring near the input end of the interferometer. The KS2 is the recommended mount for the SA200 since it will hold the SA200 tightly and provide the kinematic control necessary for proper alignment in the optical cavity. In contrast, the SA210 has a Ø1" mounting ring near the input end of the interferometer, and as a result, the recommended mount is the KS1.

SA200
Mounted in a KS2




See Page 224

MODEL	SA200	SA210
FSR ^a	1.5 GHz	10 GHz
Finesse ^b	200 (250 Typ)	150 (180 Typ)
Resolution	7.5 MHz	67 MHz
Max. Beam Diameter ^c	600 µm	150 µm
Cavity Length	50 mm	7.5 mm

^a Free spectral range for a confocal cavity and defined by $FSR = c/4r$ where c is the speed of light.

^b Effective finesse defined by $F_t = FSR/\Delta$, where $FSR = c/4r$ and $\Delta = FWHM$ impulse response.

^c Maximum beam diameter along the length of the cavity to obtain the typical measured finesse specified.



1.5 GHz Free Spectral Range, 50 mm Cavity Length, Finesse $F_t \geq 200$ (Typical $F_t = 250$)

NEW

ITEM#	\$	£	€	RMB	DESCRIPTION	RECOMMENDED MOUNT
SA200-2A	\$ 2,727.00	£ 1,890.50	€ 2,421.00	¥ 23,027.00	Scanning Fabry-Perot, 370-410 nm, 1.5 GHz FSR	KS2 (See Page 224)
SA200-5A	\$ 2,535.10	£ 1,757.50	€ 2,250.50	¥ 21,407.00	Scanning Fabry-Perot, 525-650 nm, 1.5 GHz FSR	
SA200-6A	\$ 2,828.00	£ 1,960.50	€ 2,510.50	¥ 23,880.00	Scanning Fabry-Perot, 650-800 nm, 1.5 GHz FSR	
SA200-7A	\$ 2,747.20	£ 1,904.50	€ 2,439.00	¥ 23,198.00	Scanning Fabry-Perot, 780-930 nm, 1.5 GHz FSR	
SA200-9A	\$ 2,939.10	£ 2,037.50	€ 2,609.50	¥ 24,818.00	Scanning Fabry-Perot, 900-1100 nm, 1.5 GHz FSR	
SA200-12A	\$ 2,979.50	£ 2,065.50	€ 2,645.00	¥ 25,159.00	Scanning Fabry-Perot, 1250-1400 nm, 1.5 GHz FSR	
SA200-14A	\$ 2,646.20	£ 1,834.50	€ 2,349.50	¥ 22,345.00	Scanning Fabry-Perot, 1450-1625 nm, 1.5 GHz FSR	

10 GHz Free Spectral Range, 7.5 mm Cavity Length, Finesse $F_t \geq 150$ (Typical $F_t = 180$)

ITEM#	\$	£	€	RMB	DESCRIPTION	RECOMMENDED MOUNT
SA210-5A	\$ 2,656.30	£ 1,841.50	€ 2,358.50	¥ 22,430.00	Scanning Fabry-Perot, 525-650 nm, 10 GHz FSR	KS1 (See Page 224)
SA210-6A	\$ 2,838.10	£ 1,967.50	€ 2,519.50	¥ 23,965.00	Scanning Fabry-Perot, 650-800 nm, 10 GHz FSR	
SA210-7A	\$ 2,716.90	£ 1,883.50	€ 2,412.00	¥ 22,942.00	Scanning Fabry-Perot, 780-930 nm, 10 GHz FSR	
SA210-9A	\$ 2,706.80	£ 1,876.50	€ 2,403.00	¥ 22,857.00	Scanning Fabry-Perot, 900-1100 nm, 10 GHz FSR	
SA210-12A	\$ 2,979.50	£ 2,065.50	€ 2,645.00	¥ 25,159.00	Scanning Fabry-Perot, 1250-1400 nm, 10 GHz FSR	
SA210-14A	\$ 2,747.20	£ 1,904.50	€ 2,439.00	¥ 23,198.00	Scanning Fabry-Perot, 1450-1625 nm, 10 GHz FSR	

Scanning Fabry-Perot Controller Box

Photo Amplifier Specifications

- **Gain Steps:** 0, 10, and 20 dB
- **Transimpedance Gain (Hi-Z):** 10 kV/A, 100 kV/A, and 1 MV/A
- **Transimpedance Gain (50 Ω):** 5 kV/A, 50 kV/A, and 500 kV/A
- **Output Voltage:** 0 - 10 V (Min)
- **Bandwidth:** 250 kHz
- **Noise (RMS):**
 - 0.1 mV @ 10 kV/A Gain
 - 0.2 mV @ 100 kV/A Gain
 - 1.5 mV @ 1 MV/A Gain

Ramp Specifications


- **Waveform:** Sawtooth or Triangle
- **Output Voltage Range:** 1 - 45 V (Offset + Amplitude)
- **Offset Range:** 0 - 15 VDC
- **Amplitude Range:** 1-30 V
- **Risetime Range:** 0.01 - 0.1 s (1X Sweep Expansion) 1 - 10 s (100X Sweep Expansion)
- **Sweep Expansion:** 1X, 2X, 5X, 10X, 20X, 50X, or 100X
- **Sweep Scale Error:** $\pm 0.5\%$
- **Output Noise:** 1 mV_{RMS} (~6.6 mV_{PP})
- **Trigger:** Ramp Start or Midpoint

The SA201 is specifically designed to control Thorlabs' Fabry-Perot Interferometers by generating a high-stability, low-noise voltage ramp. This ramp signal is used to scan the separation between the two cavity mirrors. The controller adjusts the ramp voltage and scan time, allowing the user to choose the scan range and speed, while an offset control allows the spectrum displayed on an oscilloscope to be shifted right or left.



A TTL output allows the user to externally trigger an oscilloscope on either the beginning or midpoint of the ramp waveform. The ability to trigger the oscilloscope from the mid-point makes zooming in on a lineshape more convenient. Simply place the spectral component of interest on the center of the screen, and increase the timebase of the oscilloscope. There is no need to use the offset to re-center the signal since the scope expands the time scale about the point of interest. The controller also has a calibrated zoom capability that provides a 1X, 2X, 5X, 10X, 20X, 50X, or 100X increase in the period of the ramp signal, which allows for an extremely wide range of scan times.


The SA201 also includes a high-precision photodetector amplifier circuit used to monitor the transmission of the cavity. The amplifier provides an adjustable transimpedance gain of 10 kV/A, 100 kV/A, or 1 MV/A when driving a high impedance load, such as an oscilloscope. Using the output sync signal from the controller, an oscilloscope can be used to display the spectrum of the input laser. The detector circuitry incorporates a blanking circuit that disables the photodiode response during the falling edge of the sawtooth waveform. The blanking circuit can be disabled by switching a circuit board jumper as described in the manual.



ITEM#	\$	£	€	RMB	DESCRIPTION
SA201	\$ 841.50	£ 583.40	€ 747,10	¥ 7,105.70	Control Box for the SA200 Family of Interferometers

Compact CCD Spectrometers

NEW
product



Features

- Wavelength Range: 3 Models Cover 200-1000 nm
- Resolution: <math><0.5\text{ nm}</math> FWHM
- Sensitivity: 160 V/lux-s
- Integration Time: 10 μs to 60 s
- Czerny-Turner Spectrometer
- Compact Size: 30 mm x 120 mm x 80 mm
- High-Speed USB Connection
- External Trigger Synchronization
- 16-Bit A/D Converter
- 3,648 Pixel CCD Line Array

See Pages 1310-1311