

# 3D Reflection from a Mirror




Expressions for the direction of the reflected ray and points on the reflected beam path.

◆ Overview .....	3
◆ Reflection from a Surface, in Local Coordinates .....	4
◆ Defining Global Coordinate System and Rotation Conventions .....	6
◆ Calculating a Reflected Ray's Global Coordinates .....	8
– Rotation Matrices .....	9
– Rotation Order Matters, Since Rotations do not Commute .....	11
– Steps of the Procedure .....	13
◆ Selected Cases of Rotating a Mounted Mirror	
– Example 1: Mount Adjusters Used to Tune both Pitch and Yaw .....	14
– Example 2: Post Rotation Used for Yaw, Mount Adjuster for Pitch .....	29

# Tracing the Reflected Beam Path

- ◆ Each mirror in a setup has its own, independently adjustable, angular orientation.
- ◆ The beam path depends on each mirror's orientation.
- ◆ Defining a fixed (global) coordinate system for the setup is useful for tracing the beam path through the setup.
- ◆ However, it is simplest to calculate the direction of the reflected beam when working in the local coordinate system of the reflective surface.
- ◆ Therefore, both local and global coordinate systems are often used, and it is necessary to convert between them.



**Figure 1.** The beam path reflected by these two mirrors depends on their orientations with respect to one another.

# Surface Reflection in Terms of Local Coordinates

The incident ray reflects across the surface normal.

- ◆ The angles of incidence and reflection with the surface normal are the same.


$$\alpha_i = \alpha_r$$

- ◆ The optical angle ( $\theta$ ) between the incident and reflected rays is twice the angle between the incident ray and surface normal.

$$\theta = 2\alpha_i$$

- ◆ Changing the angle of incidence by an angle  $\delta$  changes the optical angle ( $\theta$ ) between the incident and reflected rays by an angle  $2\delta$ .

$$\theta = 2(\alpha_i + \delta) = 2\alpha_i + 2\delta$$



**Figure 2.** Incident and reflected ray angles with the surface normal are equal.

# Surface Reflection in Terms of Local Coordinates

Calculating the reflected ray's coordinates.


- ◆ In the local coordinate system, the surface is in the plane of the  $u'$ - and  $v'$ -axes, while the  $w'$ -axis is normal to the surface and  $u'$ - $v'$  plane.
- ◆ The incident and reflected rays have unit vectors,  $i'$  and  $r'$ , respectively, in which:

$$i'_{\parallel} = r'_{\parallel} \quad \text{Components Parallel to } u'\text{-}v' \text{ Plane}$$

$$i'_{\perp} = -r'_{\perp} \quad \text{Components Perpendicular to } u'\text{-}v' \text{ Plane}$$

- ◆ The reflected ray ( $r'$ ) in local coordinates, is calculated by reflecting the incident ray ( $i'$ ) across the surface normal ( $n'$ ):

$$r' = i' - 2(i' \cdot n')n'$$




**Figure 3.** The angle between the incident ray and surface normal equals the angle between the reflected ray and surface normal.


# The Different Local and Global Perspectives

Identifying differences in and relating local and global system perspectives.

- ◆ Local perspective: the surface never moves. Instead, the angle of the incident ray changes.
- ◆ Global perspective: the surface rotates relative to the incident ray.
- ◆ Example: choose a global coordinate system and define the orientation of the unrotated reflective surface.
  - Position of surface center:  $(0, 0, 0)$ .
  - The unrotated surface is in the  $x$ - $y$  plane.
  - The  $z$ -axis is normal to the unrotated surface.



**Figure 4.** Local Perspective of Surface




**Figure 5.** First step to a global view is defining the orientation of the unrotated surface.

# Define Some Global Coordinate System Conventions

Define global coordinate conventions for rotating the surface.

- ◆ Global coordinates include information about the rotation angles of the surface relative to the global axes.
- ◆ Pitch and Yaw Rotation of the Reflective Surface
  - Positive pitch ( $\theta$ ) around x-axis is counterclockwise (CCW), when looking towards the origin, down the x-axis.
  - Positive yaw ( $\psi$ ) around y-axis is CCW, when looking towards the origin, down the y-axis.
  - The center of the reflective surface always coincides with the origin, regardless of its orientation.




**Figure 6.** Rotation angles with respect to the x- and y-axis ( $\theta$  and  $\psi$ , respectively) are measured counterclockwise.

# Reflection Matrices: Coordinates and Unit Vectors

Points on the surface have local and global coordinates and unit vectors.

- ◆ Both points and vectors can be converted between local and global coordinate systems.
- ◆ Unit vectors directed from the origin, towards a point:
  - Local unit vector:  $\mathbf{r}' = \langle u', v', w' \rangle$
  - Global unit vector:  $\mathbf{r} = \langle x, y, z \rangle$
  - For the unrotated surface,  $\mathbf{r}' = \mathbf{r}$ .
- ◆ Points on the surface:
  - Local coordinates:  $P' = (u', v', w')$
  - Global coordinates:  $P = (x, y, z)$
  - For the unrotated surface,  $P' = P$ .



**Figure 7.** Rotation around the (a) x-axis and (b) y-axis. Unrotated surface is on the left, rotated surface is on the right.

# Reflection Matrices: Unrotated to Rotated Orientations

Matrix algebra can be used to rotate vectors and points around an axis.

- ◆ Converting a local point or unit vector to global coordinates requires including information about the rotation angles. This can be done using matrices.

- ◆ When rotation is CCW around the x-axis:


$$P = R_x(\theta)P'$$

$$\begin{bmatrix} x \\ y \\ z \end{bmatrix} = \begin{bmatrix} 1 & 0 & 0 \\ 0 & \cos\theta & -\sin\theta \\ 0 & \sin\theta & \cos\theta \end{bmatrix} \begin{bmatrix} u' \\ v' \\ w' \end{bmatrix}$$

- ◆ When rotation is CCW around the y-axis:

$$P = R_y(\psi)P'$$

$$\begin{bmatrix} x \\ y \\ z \end{bmatrix} = \begin{bmatrix} \cos\psi & 0 & \sin\psi \\ 0 & 1 & 0 \\ -\sin\psi & 0 & \cos\psi \end{bmatrix} \begin{bmatrix} u' \\ v' \\ w' \end{bmatrix}$$




**Figure 8.** Counterclockwise (CCW) rotation around the (a) x-axis and (b) y-axis. Unrotated surface is on the left, rotated surface is on the right.

# Object Orientation Depends on Order of Rotations

Note that when rotating an object around a coordinate system's axes, the object's final orientation depends on the order in which the rotations were performed.


## Case 1:

### Yaw then Pitch



## Case 2:

### Pitch then Yaw



# Object Orientation Depends on Order of Rotations

In other words, rotation matrices are not commutative.

- ◆ The final orientation of the reflective surface depends on the order in which the mirror's pitch and yaw axes were adjusted.
- ◆ Therefore, the order in which the pitch and yaw axes are adjusted determines:
  - The orientation (direction) of the reflected beam.
  - The beam path.
- ◆ To ensure agreement between experimental and modeled results, the order, direction, and magnitude of the rotations in the experimental and modeled cases must be in perfect agreement.

# Reflection Matrices: Sequence of Rotations

A single, custom matrix converts between local and global coordinate systems.

- ◆ A sequence of rotations is typically used to orient a mirror. A total matrix ( $\mathbf{R}_{Total}$ ) converts between coordinate systems and accounts for all applied rotations.
- ◆ Compute the total matrix by multiplying the individual rotation matrices together. Multiply them in the order in which the sequence of rotations was performed.
- ◆ For example, if the first rotation was around the  $x$ -axis ( $\mathbf{R}_x(\theta)$ ), and the second was around the  $y$ -axis ( $\mathbf{R}_y(\psi)$ ), the rotation matrix ( $\mathbf{R}_{yx}(\theta, \psi)$ ) is the product:

$$\mathbf{P} = \mathbf{R}_{Total}\mathbf{P}' = \mathbf{R}_y(\psi)\mathbf{R}_x(\theta)\mathbf{P}' = \mathbf{R}_{yx}(\theta, \psi)\mathbf{P}'$$
$$\mathbf{R}_{yx}(\theta, \psi) = \begin{bmatrix} \cos\psi & 0 & \sin\psi \\ 0 & 1 & 0 \\ -\sin\psi & 0 & \cos\psi \end{bmatrix} \begin{bmatrix} 1 & 0 & 0 \\ 0 & \cos\theta & -\sin\theta \\ 0 & \sin\theta & \cos\theta \end{bmatrix} = \begin{bmatrix} \cos\psi & \sin\theta\sin\psi & \cos\theta\sin\psi \\ 0 & \cos\theta & -\sin\theta \\ -\sin\psi & \sin\theta\cos\psi & \cos\theta\cos\psi \end{bmatrix}$$

For information about matrix rotations and transformations, refer to a linear algebra reference, such as: D. Cherney, T. Denton, R. Thomas, and A. Waldron (Davis, CA 2013). *Linear Algebra*. <https://www.math.ucdavis.edu/~linear/linear-guest.pdf>

# Calculate the Reflected Ray's Global Coordinates

Procedure for calculating the reflected ray's direction relative to the incident ray's.

- ◆ Transform the unit vector of the incident ray from global ( $\mathbf{i}$ ) to local ( $\mathbf{i}'$ ) coordinates using the inverse rotation matrix ( $\mathbf{R}_{Total}^{-1}$ ):


$$\mathbf{i}' = \mathbf{R}_{Total}^{-1}(\theta, \psi)\mathbf{i}$$

- ◆ Reflect the ray across the local normal ( $\mathbf{n}'$ ) to obtain the reflected ray ( $\mathbf{r}'$ ) in local coordinates:

$$\mathbf{r}' = \mathbf{i}' - 2(\mathbf{i}' \cdot \mathbf{n}')\mathbf{n}'$$

- ◆ Convert the reflected ray back into global coordinates ( $\mathbf{r}$ ) using the rotation matrix ( $\mathbf{R}_{total}$ ):

$$\mathbf{r} = \mathbf{R}_{Total}(\theta, \psi)\mathbf{r}'$$



**Figure 9.** Reflection is performed using local coordinates, which involves reflecting the incident ray across the normal to the  $u$ - $v$  plane.


Example 1: Mount's Adjusters Tune the Reflected Beam's Direction

**THORLABS**

# Example 1 Overview: Adjusters Tune Orientation

Find the reflected ray's direction, relative to the incident ray, after angle tuning.

- ◆ The mount's back plate does not move.
- ◆ The mirror is installed in the mount's front plate, which can rotate around the pivot point.
- ◆ The mirror's orientation is tuned using only the mount's adjusters, which rotate the mount's front plate about the pivot point.
- ◆ The mirror rotates relative to the fixed global  $x$ -,  $y$ -, and  $z$ -axes, whose origin is chosen to be the front plate's pivot point.
- ◆ The incident ray's direction is fixed and parallel to the  $z$ -axis.




**Figure 10.** Mirror mounted in a KM200 kinematic mirror mount.

# Example 1 Overview: Adjusters Tune Orientation

Two rotations, one pitch and one yaw, were performed in succession.

- ◆ The mount's back plate cannot move, since it is secured to a post, which is clamped in a post holder.
- ◆ The pitch and yaw adjusters are installed in the back plate.
- ◆ The adjusters' tips push against the backside of the front plate.
  - Tuning the pitch adjuster rotates the front plate around the  $x$ -axis.
  - Tuning the yaw adjuster rotates the front plate around the  $y$ -axis.



**Figure 11.** Adjusters can be used to tune the mirror's pitch and yaw.

# Calculate the Reflected Ray's Global Coordinates

Procedure for calculating the reflected ray's direction relative to the incident ray's.

- ◆ Transform the unit vector of the incident ray from global ( $\mathbf{i}$ ) to local ( $\mathbf{i}'$ ) coordinates using the inverse rotation matrix ( $\mathbf{R}_{Total}^{-1}$ ):


$$\mathbf{i}' = \mathbf{R}_{Total}^{-1}(\theta, \psi)\mathbf{i}$$

- ◆ Reflect the ray across the local normal ( $\mathbf{n}'$ ) to obtain the reflected ray ( $\mathbf{r}'$ ) in local coordinates:

$$\mathbf{r}' = \mathbf{i}' - 2(\mathbf{i}' \cdot \mathbf{n}')\mathbf{n}'$$

- ◆ Convert the reflected ray back into global coordinates ( $\mathbf{r}$ ) using the rotation matrix ( $\mathbf{R}_{total}$ ):

$$\mathbf{r} = \mathbf{R}_{Total}(\theta, \psi)\mathbf{r}'$$



**Figure 12.** Reflection is performed using local coordinates, which involves reflecting the incident ray across the normal to the  $u$ - $v$  plane.

# Mount Adjusters Tuned: Local to Global Transformations

- ◆ The total rotation matrix ( $\mathbf{R}_{Total}(\theta, \psi)$ ), which converts **local** (unrotated) to **global** (rotated) coordinates, differs depending on whether pitch or yaw is adjusted first.

$$\mathbf{P}_{Global} = \mathbf{R}_{Total}(\theta, \psi) \mathbf{P}'_{Local}$$

If **yaw**, then **pitch**, was adjusted, the local to global transformation:

$$\mathbf{R}_{xy}(\theta, \psi) = \mathbf{R}_x(\theta) \mathbf{R}_y(\psi)$$

$$\begin{array}{cc} \text{Pitch} & \text{Yaw} \\ = \begin{bmatrix} 1 & 0 & 0 \\ 0 & \cos\theta & -\sin\theta \\ 0 & \sin\theta & \cos\theta \end{bmatrix} \begin{bmatrix} \cos\psi & 0 & \sin\psi \\ 0 & 1 & 0 \\ -\sin\psi & 0 & \cos\psi \end{bmatrix} \end{array}$$

$$\mathbf{R}_{xy}(\theta, \psi) = \begin{bmatrix} \cos\psi & 0 & \sin\psi \\ \sin\theta \sin\psi & \cos\theta & -\sin\theta \cos\psi \\ -\cos\theta \sin\psi & \sin\theta & \cos\theta \cos\psi \end{bmatrix}$$

If **pitch**, then **yaw**, was adjusted, the local to global transformation:

$$\mathbf{R}_{yx}(\psi, \theta) = \mathbf{R}_y(\psi) \mathbf{R}_x(\theta)$$

$$\begin{array}{cc} \text{Yaw} & \text{Pitch} \\ = \begin{bmatrix} \cos\psi & 0 & \sin\psi \\ 0 & 1 & 0 \\ -\sin\psi & 0 & \cos\psi \end{bmatrix} \begin{bmatrix} 1 & 0 & 0 \\ 0 & \cos\theta & -\sin\theta \\ 0 & \sin\theta & \cos\theta \end{bmatrix} \end{array}$$

$$\mathbf{R}_{yx}(\psi, \theta) = \begin{bmatrix} \cos\psi & \sin\theta \sin\psi & \cos\theta \sin\psi \\ 0 & \cos\theta & -\sin\theta \\ -\sin\psi & \sin\theta \cos\psi & \cos\theta \cos\psi \end{bmatrix}$$

# Mount Adjusters Tuned: Global to Local Transformations

The inverse total rotation matrix ( $\mathbf{R}_{Total}^{-1}(\theta, \psi)$ ) converts global to local coordinates.

- ◆ **Global** (rotated) to **local** (unrotated) coordinates:  $\mathbf{P}'_{Local} = \mathbf{R}^{-1}(\theta, \psi)\mathbf{P}_{Global}$ 
  - $\mathbf{R}^{-1}(\theta, \psi)$  is the inverse of  $\mathbf{R}(\theta, \psi)$ .
  - In the case of rotation matrices, the inverse equals the transpose.

If **yaw**, then **pitch**, was adjusted, the global to local transformation:

$$\mathbf{R}_{xy}^{-1}(\theta, \psi) = \begin{bmatrix} \cos \psi & \sin \theta \sin \psi & -\cos \theta \sin \psi \\ 0 & \cos \theta & \sin \theta \\ \sin \psi & -\sin \theta \cos \psi & \cos \theta \cos \psi \end{bmatrix}$$

If **pitch**, then **yaw**, was adjusted, the global to local transformation:

$$\mathbf{R}_{yx}^{-1}(\psi, \theta) = \begin{bmatrix} \cos \psi & 0 & -\sin \psi \\ \sin \theta \sin \psi & \cos \theta & \sin \theta \cos \psi \\ \cos \theta \sin \psi & -\sin \theta & \cos \theta \cos \psi \end{bmatrix}$$

# Mount Adjusters Tuned: Incident Ray in Local Coords

Start by transforming the unit vector of the incident ray into local coordinates.

- ◆ Since incident ray is parallel to the z-axis, the unit vector of the incident ray in global coordinates is:  $\mathbf{i} = \langle 0, 0, -1 \rangle = \begin{bmatrix} 0 \\ 0 \\ -1 \end{bmatrix}$

If **yaw**, then **pitch**, was adjusted, the incident ray in local coordinates:

$$\mathbf{i}'_{Local} = \mathbf{R}_{xy}^{-1}(\theta, \psi) \mathbf{i}_{Global}$$

$$\begin{bmatrix} u' \\ v' \\ w' \end{bmatrix} = \begin{bmatrix} \cos \psi & \sin \theta \sin \psi & -\cos \theta \sin \psi \\ 0 & \cos \theta & \sin \theta \\ \sin \psi & -\sin \theta \cos \psi & \cos \theta \cos \psi \end{bmatrix} \begin{bmatrix} 0 \\ 0 \\ -1 \end{bmatrix}$$

$$\begin{bmatrix} u' \\ v' \\ w' \end{bmatrix} = \begin{bmatrix} +\cos \theta \sin \psi \\ -\sin \theta \\ -\cos \theta \cos \psi \end{bmatrix}$$

If **pitch**, then **yaw**, was adjusted, the incident ray in local coordinates:

$$\mathbf{i}'_{Local} = \mathbf{R}_{yx}^{-1}(\psi, \theta) \mathbf{i}_{Global}$$

$$\begin{bmatrix} u' \\ v' \\ w' \end{bmatrix} = \begin{bmatrix} \cos \psi & 0 & -\sin \psi \\ \sin \theta \sin \psi & \cos \theta & \sin \theta \cos \psi \\ \cos \theta \sin \psi & -\sin \theta & \cos \theta \cos \psi \end{bmatrix} \begin{bmatrix} 0 \\ 0 \\ -1 \end{bmatrix}$$

$$\begin{bmatrix} u' \\ v' \\ w' \end{bmatrix} = \begin{bmatrix} +\sin \psi \\ -\sin \theta \cos \psi \\ -\cos \theta \cos \psi \end{bmatrix}$$

# Mount Adjusters Tuned: Reflected Ray in Local Coords

Calculate the unit vector of the reflected ray in local coordinates.

- ◆ Use the incident ray, in local coordinates, to calculate the reflected ray in local coordinates:

$$\mathbf{r}'_{Local} = \mathbf{i}'_{Local} - 2(\mathbf{i}'_{Local} \cdot \mathbf{n}'_{Local})\mathbf{n}'_{Local}$$

If **yaw**, then **pitch**, was adjusted, the reflected ray in local coordinates:

$$\begin{bmatrix} u' \\ v' \\ w' \end{bmatrix} = \begin{bmatrix} + \cos \theta \sin \psi \\ - \sin \theta \\ - \cos \theta \cos \psi \end{bmatrix} - 2 \left\{ \begin{bmatrix} + \cos \theta \sin \psi \\ - \sin \theta \\ - \cos \theta \cos \psi \end{bmatrix} \cdot \begin{bmatrix} 0 \\ 0 \\ 1 \end{bmatrix} \right\} \begin{bmatrix} 0 \\ 0 \\ 1 \end{bmatrix}$$

$$\begin{bmatrix} u' \\ v' \\ w' \end{bmatrix} = \begin{bmatrix} + \cos \theta \sin \psi \\ - \sin \theta \\ + \cos \theta \cos \psi \end{bmatrix}$$

If **pitch**, then **yaw**, was adjusted, the reflected ray in local coordinates:

$$\begin{bmatrix} u' \\ v' \\ w' \end{bmatrix} = \begin{bmatrix} + \sin \psi \\ - \sin \theta \cos \psi \\ - \cos \theta \cos \psi \end{bmatrix} - 2 \left\{ \begin{bmatrix} + \sin \psi \\ - \sin \theta \cos \psi \\ - \cos \theta \cos \psi \end{bmatrix} \cdot \begin{bmatrix} 0 \\ 0 \\ 1 \end{bmatrix} \right\} \begin{bmatrix} 0 \\ 0 \\ 1 \end{bmatrix}$$

$$\begin{bmatrix} u' \\ v' \\ w' \end{bmatrix} = \begin{bmatrix} + \sin \psi \\ - \sin \theta \cos \psi \\ + \cos \theta \cos \psi \end{bmatrix}$$

# Mount Adjusters Tuned: Reflected Ray in Global Coords

Calculate the unit vector of the reflected ray in global coordinates.

- ◆ Transform the reflected ray from local coordinates into global coordinates:

If **yaw**, then **pitch**, was adjusted,  
the reflected ray in global coordinates:

$$\mathbf{r}_{Global} = \mathbf{R}_{xy}(\theta, \psi) \mathbf{r}'_{Local}$$

$$\begin{bmatrix} x \\ y \\ z \end{bmatrix} = \begin{bmatrix} \cos \psi & 0 & \sin \psi \\ \sin \theta \sin \psi & \cos \theta & -\sin \theta \cos \psi \\ -\cos \theta \sin \psi & \sin \theta & \cos \theta \cos \psi \end{bmatrix} \begin{bmatrix} +\cos \theta \sin \psi \\ -\sin \theta \\ +\cos \theta \cos \psi \end{bmatrix}$$

$$= \begin{bmatrix} 2 \cos \theta \cos \psi \sin \psi \\ -2 \cos \theta \sin \theta \cos^2 \psi \\ -\cos^2 \theta \sin^2 \psi - \sin^2 \theta + \cos^2 \theta \cos^2 \psi \end{bmatrix}$$

If **pitch**, then **yaw**, was adjusted,  
the reflected ray in global coordinates:

$$\mathbf{r}_{Global} = \mathbf{R}_{xy}(\theta, \psi) \mathbf{r}'_{Local}$$

$$\begin{bmatrix} x \\ y \\ z \end{bmatrix} = \begin{bmatrix} \cos \psi & \sin \theta \sin \psi & \cos \theta \sin \psi \\ 0 & \cos \theta & -\sin \theta \\ -\sin \psi & \sin \theta \cos \psi & \cos \theta \cos \psi \end{bmatrix} \begin{bmatrix} +\sin \psi \\ -\sin \theta \cos \psi \\ +\cos \theta \cos \psi \end{bmatrix}$$

$$= \begin{bmatrix} \cos \psi \sin \psi (1 - \sin^2 \theta + \cos^2 \theta) \\ -2 \cos \theta \sin \theta \cos \psi \\ -\sin^2 \psi - \sin^2 \theta \cos^2 \psi + \cos^2 \theta \cos^2 \psi \end{bmatrix}$$

# Mount Adjusters Tuned: Reflected Ray in Global Coords

Simplify the z-component for the reflected ray in global coordinates.

- ◆ If **yaw**, then **pitch**, was adjusted, the reflected vector in global coordinates:

$$\begin{bmatrix} x \\ y \\ z \end{bmatrix} = \begin{bmatrix} 2 \cos \theta \cos \psi \sin \psi \\ -2 \cos \theta \sin \theta \cos^2 \psi \\ \boxed{-\cos^2 \theta \sin^2 \psi - \sin^2 \theta + \cos^2 \theta \cos^2 \psi} \end{bmatrix}$$

$$\begin{aligned} &= -\cos^2 \theta \sin^2 \psi - \sin^2 \theta + \cos^2 \theta \cos^2 \psi + (\cos^2 \theta \cos^2 \psi - \cos^2 \theta \cos^2 \psi) \\ &= -(\cos^2 \theta \{\cos^2 \psi + \sin^2 \psi\} + \sin^2 \theta) + 2 \cos^2 \theta \cos^2 \psi \\ &= -(\cos^2 \theta \{1\} + \sin^2 \theta) + 2 \cos^2 \theta \cos^2 \psi \\ &= -1 + 2 \cos^2 \theta \cos^2 \psi \end{aligned}$$

# Mount Adjusters Tuned: Reflected Ray in Global Coords

Simplify the x- and z-components for the reflected ray in global coordinates.

- ◆ If pitch, then yaw, was adjusted, the reflected vector in global coordinates:

$$\begin{bmatrix} x \\ y \\ z \end{bmatrix} = \begin{bmatrix} \cos \psi \sin \psi (1 - \sin^2 \theta + \cos^2 \theta) \\ -2 \cos \theta \sin \theta \cos^2 \psi \\ -\sin^2 \psi - \sin^2 \theta \cos^2 \psi + \cos^2 \theta \cos^2 \psi \end{bmatrix}$$

$\begin{aligned} &= \cos \psi \sin \psi (\cos^2 \theta + \cos^2 \theta) \\ &= 2 \cos^2 \theta \cos \psi \sin \psi \end{aligned}$

$\begin{aligned} &= -1 + \cos^2 \psi - \sin^2 \theta \cos^2 \psi + \cos^2 \theta \cos^2 \psi \\ &= -1 + \cos^2 \psi (\cos^2 \theta + \sin^2 \theta) - \sin^2 \theta \cos^2 \psi + \cos^2 \theta \cos^2 \psi \\ &= -1 + 2\cos^2 \theta \cos^2 \psi \end{aligned}$

# Mount Adjusters Tuned: Reflected Ray in Global Coords

Expression for the unit vector of the reflected ray in global coordinates.

- ◆ The z-component is the same in both cases, but the x- and y-components differ.

If first the mirror's **yaw ( $\psi$ )** and then its **pitch ( $\theta$ )** was adjusted, the reflected ray in global coordinates:

$$\begin{bmatrix} x \\ y \\ z \end{bmatrix} = \begin{bmatrix} 2 \cos \theta \cos \psi \sin \psi \\ -2 \cos \theta \sin \theta \cos^2 \psi \\ -1 + 2 \cos^2 \theta \cos^2 \psi \end{bmatrix}$$

If first the mirror's **pitch ( $\theta$ )** and then its **yaw ( $\psi$ )** was adjusted, the reflected ray in global coordinates:

$$\begin{bmatrix} x \\ y \\ z \end{bmatrix} = \begin{bmatrix} 2 \cos^2 \theta \cos \psi \sin \psi \\ -2 \cos \theta \sin \theta \cos \psi \\ -1 + 2 \cos^2 \theta \cos^2 \psi \end{bmatrix}$$

# Mount Adjusters Tuned: Using the Reflected Ray

Use the reflected ray's unit vector to calculate arbitrary points on the ray path.

- ◆ The reflected ray's unit vector  $\langle x, y, z \rangle$  is useful because,
  - It points in the direction of the ray's path.
  - It has a length of 1, so the point  $(x, y, z)$  lies at the ray's tip:  $\sqrt{x^2 + y^2 + z^2} = 1$ .
  - It can be used to calculate the coordinates of any point on the beam path.
- ◆ To calculate an arbitrary point  $(x_2, y_2, z_2)$  on the reflected ray's path,
  - Calculate a new vector by multiplying the unit vector by a constant  $A$ :
$$A\langle x, y, z \rangle = \langle Ax, Ay, Az \rangle$$
  - The length of the new vector is  $A$ :  $\sqrt{|Ax|^2 + |Ay|^2 + |Az|^2} = A$
  - Choose the length so that:  $(x_2, y_2, z_2) = (Ax, Ay, Az)$
  - The point  $(x_2, y_2, z_2)$  lies at the new vector's tip.

# Mount Adjusters Tuned: Points on the Reflected Ray

Calculating an arbitrary point on the beam path when another point is known.

- ◆ The known point is:  $(x_1, y_1, z_1)$ .
- ◆ The coordinates of the known point can be added to a vector equal to the reflected unit vector whose length has been scaled by the correct value of A.

If **yaw ( $\psi$ )**, then **pitch ( $\theta$ )**, was adjusted, the unit vector of the reflected ray:

$$\begin{bmatrix} x \\ y \\ z \end{bmatrix} = \begin{bmatrix} 2 \cos \theta \cos \psi \sin \psi \\ -2 \cos \theta \sin \theta \cos^2 \psi \\ -1 + 2 \cos^2 \theta \cos^2 \psi \end{bmatrix}$$

Points on the reflected ray:

$$\begin{aligned} x_2 &= x_1 + A(2 \cos \theta \cos \psi \sin \psi) \\ y_2 &= y_1 + A(-2 \cos \theta \sin \theta \cos^2 \psi) \\ z_2 &= z_1 + A(-1 + 2 \cos^2 \theta \cos^2 \psi) \end{aligned}$$

If **pitch ( $\theta$ )**, then **yaw ( $\psi$ )**, was adjusted, the unit vector of the reflected ray:

$$\begin{bmatrix} x \\ y \\ z \end{bmatrix} = \begin{bmatrix} 2 \cos^2 \theta \cos \psi \sin \psi \\ -2 \cos \theta \sin \theta \cos \psi \\ -1 + 2 \cos^2 \theta \cos^2 \psi \end{bmatrix}$$

Points on the reflected ray:

$$\begin{aligned} x_2 &= x_1 + A(2 \cos^2 \theta \cos \psi \sin \psi) \\ y_2 &= y_1 + A(-2 \cos \theta \sin \theta \cos \psi) \\ z_2 &= z_1 + A(-1 + 2 \cos^2 \theta \cos^2 \psi) \end{aligned}$$

# Mount Adjusters Tuned: Points on the Reflected Ray

Assume the known point is at the origin:  $(x_1, y_1, z_1) = (0, 0, 0)$ .

- ◆ Calculate the required spacing between paired steering mirrors.
  - Option 1: Vary the scaling factor ( $A$ ) to change the distance from origin to beam point.
  - Option 2: Use a known point coordinate (e.g. a new beam height  $y_2$ ), the pitch angle, and the yaw angle to calculate  $A$ . Then, calculate the other two point coordinates.

If **yaw ( $\psi$ )**, then **pitch ( $\theta$ )**, was adjusted, the separation between the origin and the point  $(x_2, y_2, z_2)$  on the reflected ray:

$$\begin{aligned}x_2 &= A(2 \cos \theta \cos \psi \sin \psi) \\y_2 &= A(-2 \cos \theta \sin \theta \cos^2 \psi) \\z_2 &= A(-1 + 2 \cos^2 \theta \cos^2 \psi)\end{aligned}$$

If **pitch ( $\theta$ )**, then **yaw ( $\psi$ )**, was adjusted, the separation between the origin and the point  $(x_2, y_2, z_2)$  on the reflected ray:

$$\begin{aligned}x_2 &= A(2 \cos^2 \theta \cos \psi \sin \psi) \\y_2 &= A(-2 \cos \theta \sin \theta \cos \psi) \\z_2 &= A(-1 + 2 \cos^2 \theta \cos^2 \psi)\end{aligned}$$


Example 2: Mount Rotation for Yaw, Adjuster Tuning for Pitch

**THORLABS**

## Example 2 Overview: Choose Global Coordinates

Define a fixed, global coordinate system ( $x$ -,  $y$ -, and  $z$ -axes).

- ◆ The chosen fixed and global coordinate system:
  - The  $x$ - $y$  plane is in the plane of the unrotated mirror.
  - The  $z$ -axis is aligned with the incident ray, whose orientation is fixed.
  - The global coordinate system's origin is placed at the center of the unrotated mirror's surface.
- ◆ The directions of the incident and reflected rays are defined relative to these fixed, global coordinate axes.
- ◆ Rotating the mirror moves it relative to the global coordinate system.




**Figure 13.** The  $x$ -,  $y$ -, and  $z$ -axes of the global coordinate system.

# Example 2 Overview: Rotate Entire Mount to Change Yaw

Effect of rotating the mount around the post axis.

- ◆ Rotation around the post axis is an alternative to using the yaw adjuster.
  - Both front and back plates rotate together as one rigid unit.
  - The yaw adjuster would have rotated the front plate rotate to the back plate.
- ◆ The effect of rotating the mount:
  - The reflective surface of the rotated mirror is at an angle to the x-y plane.
  - The mirror's surface is no longer normal to the incident ray and z-axis.
  - From the mount's point of view, the direction of the incident ray has changed.




**Figure 14.** Entire mount is rotated around the post axis to change yaw.

# Example 1 Overview: Use Mount's Adjuster to Tune Pitch

Procedure for and effect of changing the mirror's pitch.

- ◆ Pitch is adjusted by tuning the mount's pitch adjuster.
  - Mount's pitch adjuster is anchored in the back plate, which does not move.
  - Adjuster's tip presses against the front plate's backside, forcing it to rotate.
- ◆ Tuning pitch rotates the front plate around the  $x'$ -,  $y'$ -,  $z'$ -axes.
  - These axes are anchored to the mount's back plate, and the origin is the pivot point.
  - Front plate rotates around the axes' origin, relative to the back plate.



**Figure 15.** Mirror mounted in a KM200 kinematic mirror mount.

# Problem Statement: Combining Rotation Types

Combining mount rotation around post axis with adjuster tuning of front plate.

- ◆ Effect of rotating the mount around the post axis to provide yaw:

- Front plate's position relative to back plate remains the same.
- Incident ray's angle of incidence with the mount changes.

- ◆ Effect of using mount adjusters.

- Front plate moves relative to the back plate.
- Incident ray's angle of incidence with respect to the mount does not change.

- ◆ Reflected ray's final direction is the same, whether post axis rotation precedes or follows adjuster tuning.

- ◆ Include the effect of the post axis rotation by converting the incident ray's unit vector into the mount's  $(x', y', z')$  coordinates.



Figure 16. Rotation Axes

# Calculate the Reflected Ray's Global Coordinates

Convert the incident ray to mount coordinates to account for post axis rotation.

- ◆ Convert incident ray's unit vector from global ( $\mathbf{i}$ ) to mount ( $\mathbf{i}'$ ) coordinates using the ( $\mathbf{R}_y^{-1}$ ) matrix:

$$\mathbf{i}' = \mathbf{R}_y^{-1}(\phi)\mathbf{i}$$

- ◆ Convert incident ray from mount ( $\mathbf{i}'$ ) to mirror ( $\mathbf{i}''$ ) coordinates using inverse rotation matrix ( $\mathbf{R}_{Total}^{-1}$ ):

$$\mathbf{i}'' = \mathbf{R}_{Total}^{-1}(\theta, \psi)\mathbf{i}'$$

- ◆ Reflect the ray across the local normal ( $\mathbf{n}''$ ) to obtain the reflected ray ( $\mathbf{r}''$ ) in local coordinates:

$$\mathbf{r}'' = \mathbf{i}'' - 2(\mathbf{i}'' \cdot \mathbf{n}'')\mathbf{n}''$$

- ◆ Convert reflected ray from mirror ( $\mathbf{r}''$ ) to mount ( $\mathbf{r}'$ ) coordinates using the rotation matrix ( $\mathbf{R}_{Total}$ ):

$$\mathbf{r}' = \mathbf{R}_{Total}(\theta, \psi)\mathbf{r}''$$

- ◆ Convert the reflected ray from mount ( $\mathbf{r}'$ ) to global ( $\mathbf{r}$ ) coordinates using the rotation matrix ( $\mathbf{R}_y$ ):

$$\mathbf{r} = \mathbf{R}_y(\phi)\mathbf{r}'$$

# Post Axis Rotation Combined with Pitch Adjuster Tuning

Example: Reflected ray orientated via post axis rotation and mount adjuster tuning.

- ◆ Global coordinates of incident ray's unit vector:  $\mathbf{i} = \langle 0, 0, -1 \rangle = \begin{bmatrix} 0 \\ 0 \\ -1 \end{bmatrix}$
- ◆ Express the orientation of the incident ray's unit vector in the mount's coordinate system. Use the inverse yaw rotation matrix ( $\mathbf{R}_y^{-1}(\phi)$ ), where  $\phi$  is the CCW rotation angle of the mount around the post axis.

$$\mathbf{i}'_{Mount} = \mathbf{R}_y^{-1}(\phi) \mathbf{i}_{Global}$$

$$\begin{bmatrix} x' \\ y' \\ z' \end{bmatrix} = \begin{bmatrix} \cos \phi & 0 & -\sin \phi \\ 0 & 1 & 0 \\ \sin \phi & 0 & \cos \phi \end{bmatrix} \begin{bmatrix} 0 \\ 0 \\ -1 \end{bmatrix} = \begin{bmatrix} \sin \phi \\ 0 \\ -\cos \phi \end{bmatrix}$$

# Post Axis Rotation Combined with Pitch Adjuster Tuning

Example: Reflected ray orientated via post axis rotation and mount adjuster tuning.

- ◆ Express the incident ray in the local coordinates of the reflective surface, as was done in Example 1.
- ◆ In this case, only the pitch adjuster was tuned ( $\mathbf{R}_{Total}^{-1}(\theta, \psi) = \mathbf{R}_x^{-1}(\theta)$ ). Tuning the adjuster rotates the mirror CCW around the  $x'$ -axis by an angle  $\theta$ .

$$\mathbf{i}''_{Local} = \mathbf{R}_{Total}^{-1}(\theta, \psi) \mathbf{i}'_{Mount} = \mathbf{R}_x^{-1}(\theta) \mathbf{i}'_{Mount}$$

$$\begin{bmatrix} u' \\ v' \\ w' \end{bmatrix} = \begin{bmatrix} 1 & 0 & 0 \\ 0 & \cos \theta & \sin \theta \\ 0 & -\sin \theta & \cos \theta \end{bmatrix} \begin{bmatrix} \sin \phi \\ 0 \\ -\cos \phi \end{bmatrix} = \begin{bmatrix} \sin \phi \\ -\sin \theta \cos \phi \\ -\cos \theta \cos \phi \end{bmatrix}$$

# Post Axis Rotation Combined with Pitch Adjuster Tuning

Example: Reflected ray orientated via post axis rotation and mount adjuster tuning.

- ◆ Reflect the incident ray across the surface normal. The result is the unit vector of the reflected ray, expressed in the local coordinates of the mirror.

$$\mathbf{r}'' = \mathbf{i}'' - 2(\mathbf{i}'' \cdot \mathbf{n}'')\mathbf{n}''$$

$$\begin{bmatrix} u' \\ v' \\ w' \end{bmatrix} = \begin{bmatrix} + \sin \phi \\ - \sin \theta \cos \phi \\ - \cos \theta \cos \phi \end{bmatrix} - 2 \left\{ \begin{bmatrix} + \sin \phi \\ - \sin \theta \cos \phi \\ - \cos \theta \cos \phi \end{bmatrix} \cdot \begin{bmatrix} 0 \\ 0 \\ 1 \end{bmatrix} \right\} \begin{bmatrix} 0 \\ 0 \\ 1 \end{bmatrix}$$

$$\begin{bmatrix} u' \\ v' \\ w' \end{bmatrix} = \begin{bmatrix} + \sin \phi \\ - \sin \theta \cos \phi \\ + \cos \theta \cos \phi \end{bmatrix}$$

# Post Axis Rotation Combined with Pitch Adjuster Tuning

Example: Reflected ray orientated via post axis rotation and mount adjuster tuning.

- ◆ Express the reflected ray in the coordinates of the mount. Since only the pitch adjuster was tuned, ( $\mathbf{R}_{Total}(\theta, \psi) = \mathbf{R}_x(\theta)$ ).

$$\mathbf{r}' = \mathbf{R}_{total}(\theta, \psi)\mathbf{r}'' = \mathbf{R}_x(\theta)\mathbf{r}''$$

$$\begin{bmatrix} x' \\ y' \\ z' \end{bmatrix} = \begin{bmatrix} 1 & 0 & 0 \\ 0 & \cos \theta & -\sin \theta \\ 0 & \sin \theta & \cos \theta \end{bmatrix} \begin{bmatrix} + \sin \phi \\ - \sin \theta \cos \phi \\ + \cos \theta \cos \phi \end{bmatrix}$$

$$\begin{bmatrix} x' \\ y' \\ z' \end{bmatrix} = \begin{bmatrix} + \sin \phi \\ -2 \cos \theta \sin \theta \cos \phi \\ +(\cos^2 \theta - \sin^2 \theta) \cos \phi \end{bmatrix}$$

# Post Axis Rotation Combined with Pitch Adjuster Tuning

Example: Reflected ray orientated via post axis rotation and mount adjuster tuning.

- ◆ Express the unit vector of the reflected ray in global coordinates using the  $R_y(\phi)$  matrix. Compare this result to those on Slide 25.

$$\mathbf{r} = R_y(\phi)\mathbf{r}'$$

$$\begin{bmatrix} x \\ y \\ z \end{bmatrix} = \begin{bmatrix} \cos \phi & 0 & \sin \phi \\ 0 & 1 & 0 \\ -\sin \phi & 0 & \cos \phi \end{bmatrix} \begin{bmatrix} +\sin \phi \\ -2 \cos \theta \sin \theta \cos \phi \\ +(\cos^2 \theta - \sin^2 \theta) \cos \phi \end{bmatrix} = \begin{bmatrix} +\cos \phi \sin \phi + (\cos^2 \theta - \sin^2 \theta) \cos \phi \sin \phi \\ -2 \cos \theta \sin \theta \cos \phi \\ -\sin^2 \phi + (\cos^2 \theta - \sin^2 \theta) \cos^2 \phi \end{bmatrix}$$

$$\begin{bmatrix} x \\ y \\ z \end{bmatrix} = \begin{bmatrix} (\cos^2 \theta + \sin^2 \theta + \cos^2 \theta - \sin^2 \theta) \cos \phi \sin \phi \\ -2 \cos \theta \sin \theta \cos \phi \\ -1 + (\cos^2 \theta + \sin^2 \theta) \cos^2 \phi + (\cos^2 \theta - \sin^2 \theta) \cos^2 \phi \end{bmatrix}$$

$$\begin{bmatrix} x \\ y \\ z \end{bmatrix} = \begin{bmatrix} 2 \cos^2 \theta \cos \phi \sin \phi \\ -2 \cos \theta \sin \theta \cos \phi \\ -1 + 2 \cos^2 \theta \cos^2 \phi \end{bmatrix}$$

# Post Axis Rotation Combined with Pitch Adjuster Tuning

Calculate arbitrary points  $(x_2, y_2, z_2)$  on the reflected ray.

- ◆ Start with the reflected ray's unit vector:

$$\begin{bmatrix} x \\ y \\ z \end{bmatrix} = \begin{bmatrix} 2 \cos^2 \theta \cos \phi \sin \phi \\ -2 \cos \theta \sin \theta \cos \phi \\ -1 + 2 \cos^2 \theta \cos^2 \phi \end{bmatrix}$$

- ◆ Multiply the reflected ray's unit vector,  $\langle x, y, z \rangle$ , by a scaling factor ( $A$ ). Add the result to the coordinates of a known point  $(x_1, y_1, z_1)$  on the reflected ray.

$$\begin{aligned} x_2 &= x_1 + A(2 \cos^2 \theta \cos \phi \sin \phi) \\ y_2 &= y_1 + A(-2 \cos \theta \sin \theta \cos \phi) \\ z_2 &= z_1 + A(-1 + 2 \cos^2 \theta \cos^2 \phi) \end{aligned}$$

- ◆ If the reflected ray's path began at the origin  $((x_1, y_1, z_1) = (0, 0, 0))$ , points on the reflected ray can be calculated using the equations at the right and varying the scaling factor ( $A$ ).

$$\begin{aligned} x_2 &= A(2 \cos^2 \theta \cos \phi \sin \phi) \\ y_2 &= A(-2 \cos \theta \sin \theta \cos \phi) \\ z_2 &= A(-1 + 2 \cos^2 \theta \cos^2 \phi) \end{aligned}$$