

Quick-Release 60 mm Cage Mount



QRC2A
Two SM2RR Retaining
Rings Included

- Internal SM2 (Ø2.035"-40) Threading
- Dovetail Design Provides Smooth Release and Clamping Action
- One-Handed Operation for Removal and Repositioning of an Optic
- Accepts Ø2" (Ø50.8 mm) Optics up to 0.15" (3.8 mm) Thick

The QRC2A cage mount is equipped with a spring-loaded dovetail clamp that can be used to secure or release the cage mount from the two upper ER rods that it attaches to in the cage structure. The cage plate has internal SM2 threading that allows for components up to 0.15" (3.8 mm) thick to be secured with the included SM2RR retaining rings.

ITEM#	\$	£	€	RMB	DESCRIPTION
QRC2A	\$ 70.00	£ 48.60	€ 62.20	¥ 591.10	SM2 (Ø2.035"-40) Series Quick-Release Cage Mount
SM2RR	\$ 6.75	£ 4.70	€ 6.00	¥ 57.00	SM2 (Ø2.035"-40) Retaining Ring

60 mm Cage Plate with Quick-Release Magnetic Plate



LCP90F-F

LCP90F-B

- Unique Dovetail Design Provides Smooth Release and Clamping Action
- One-Handed Operation for Removal and Repositioning of an Optic
- Accepts Ø2" (Ø50.8 mm) Optics up to 0.35" (8.7 mm) Thick
- SM2 (Ø2.035"-40) Internal Thread, 0.50" (12.7 mm) Thick

The LCP90F is a magnetically coupled 60 mm Cage Plate set. The three-point kinematic design offers excellent stability, precision, and repeatability. Once the back plate is locked in position, the removable front plate can be easily inserted or removed from a fully assembled cage system by rotating it sideways and then extracting it from the cage. The internal SM2 (2.035"-40) threading allows components up to 0.35" (8.7 mm) thick to be secured into the front plate with the included SM2RR retaining ring. In addition, the 0.31" (7.9 mm) thick back plate contains SM2 threads ideal for attaching SM2 Series Lens Tubes, components, and adapters. Additional front (LCP90F-F) and back (LCP90F-B) plates are available to easily swap filters or optics.



ITEM#	\$	£	€	RMB	DESCRIPTION
LCP90F	\$ 58.00	£ 40.30	€ 51.50	¥ 489.80	60 mm Removable Cage Plate, Front and Back Plates
LCP90F-B	\$ 32.00	£ 22.20	€ 28.50	¥ 270.30	60 mm Removable Cage Plate, Back Plate
LCP90F-F	\$ 32.00	£ 22.20	€ 28.50	¥ 270.30	60 mm Removable Cage Plate, Front Plate

Nikon F-Mount Adapter for 60 mm Cage Systems



LCP04

One SM2RR Retaining Ring Included

- Attaches Nikon F-Mount to 60 mm Cage Assemblies
- #8-32 (M4 x 0.7) Tapped Holes for Post Mounting
- Spring-Actuated Lock/Release Button
- SM2 (Ø2.035"-40) Internal Thread

See Page 309 for Other Camera Mounts

The LCP04 is a 60 mm cage plate F-Mount adapter. It has internal SM2 threading on the back side of the adapter for lens tube integration. For post mounting, there are #8-32 and M4 x 0.7 mounting holes on the side of this plate. The mount can also be used to hold optics up to 0.33" (8.4 mm) using an SM2RR Retaining Ring (included). Additional rings are available separately.



LCP04
All Components
Sold Separately

ITEM#	\$	£	€	RMB	DESCRIPTION
LCP04	\$ 133.00	£ 92.20	€ 118.10	¥ 1,123.10	Nikon F-Mount Adapter for 60 mm Cage System
SM2RR	\$ 6.75	£ 4.70	€ 6.00	¥ 57.00	SM2 (Ø2.035"-40) Retaining Ring